

Physiology Lessons  
for use with the  
Biopac Student Lab

PC under Windows® 98SE, Me, 2000 Pro  
or Macintosh® 8.6 – 9.1

Manual Revision  
PL3.6.7-ML3.0.7/061903

Richard Pflanzner, Ph.D.  
*Associate Professor*  
Indiana University School of Medicine  
Purdue University School of Science

J.C. Uyehara, Ph.D.  
*Biologist*  
BIOPAC Systems, Inc.

William McMullen  
*Vice President*  
BIOPAC Systems, Inc.

**BIOPAC Systems, Inc.**



42 Aero Camino, Goleta, CA 93117  
(805) 685-0066, Fax (805) 685-0067

Email: [info@biopac.com](mailto:info@biopac.com)

Web Site: <http://www.biopac.com>



**Lesson 1 Data Report**  
**ELECTROMYOGRAPHY I**  
*Standard and Integrated EMG*



# ELECTROMYOGRAPHY I

## Standard and Integrated EMG

### DATA REPORT

Student's Name: 解正平、劉維凱、張景程

Lab Section: MD303 Group 3

Date: 2018/3/23

#### I. Data and Calculations

##### Subject Profile

Name 張景程

Height 178cm

Age 21

Weight 67kg

Gender: ☒ Male / Female

##### A. EMG Measurements

Cluster #	Forearm 1 (Dominant)				Forearm 2			
	Min [3 min]	Max [3 max]	P-P [3 p-p]	Mean [40 mean]	Min [3 min]	Max [3 max]	P-P [3 p-p]	Mean [40 mean]
1	-0.13977	0.10925	0.24902	0.02701	-0.14282	0.10742	0.25024	0.02504
2	-0.19403	0.19470	0.38513	0.05098	-0.18433	0.14587	0.33020	0.03544
3	-0.40466	0.34119	0.74585	0.07619	-0.20569	0.21362	0.41931	0.05761
4	-0.48462	0.51514	0.99976	0.11153	-0.37109	0.42175	0.79285	0.08358

Note: "Clusters" are the EMG bursts associated with each clench.

- B. Use the mean measurement from the table above to compute the percentage increase in EMG activity recorded between the weakest clench and the strongest clench of Forearm 1.

Calculation:

$$(0.11153 - 0.02701) / 0.02701 = 3.1292$$

Answer: 312.92 %

## C. Tonus Measurements

Cluster #	<i>Forearm 1 (Dominant)</i>		<i>Forearm 2</i>	
	P-P [3 p-p]	Mean [40 mean]	P-P [3 p-p]	Mean [40 mean]
1	0.04028	0.00890	0.00916	0.00146
2	0.03601	0.00411	0.01526	0.00196
3	0.04028	0.00717	0.02319	0.00314
4	0.02441	0.00741	0.03418	0.00517

## II. QUESTIONS

- D. Compare the mean measurement for the right and left maximum clench EMG cluster. Are they the same or different?

\_\_\_\_ Same ☒ Different

Which one suggests the greater clench strength?

☒ Right \_\_\_\_ Left \_\_\_\_ Neither

Explain.

一般而言，在日常生活中慣用手的使用率和鍛煉會比非慣用手要來得多，以至於慣用手的肌肉會較為發達、肌纖維數量較多，且肌肉的運動是由於肌肉纖維群受到刺激，以電位差的形式傳遞訊號，因此具有較多肌纖維的慣用手離子交換量越大，所產生的總電位差也較高，也能產生較大的力量，和我們實驗結果也較為符合。

- E. What factors in addition to sex contribute to observed differences in clench strength?

1. 年齡：尚未發育完全的兒童和青少年，其肌肉纖維較少，所以肌力也會比較少；此外，當人從中壯年逐漸步入老年時，肌力也會隨著年齡的增加而逐漸下降。 2. 肌肉組成：包含肌纖維的面積、肌肉纖維收縮時的長度、肌肉收縮時所牽扯的肌纖維數量、肌腱連接的方式等等也會影響肌力的表現。 3. 受試者的心理狀態：如果受測者處在壓力、疲勞、緊張、肌肉傷害的狀態下，也會影響測試結果。 4. 其他：環境的溫度、濕度、體溫、皮膚含水量等等因素也會影響測量時所量測到的肌肉電位。

- F. Does there appear to be any difference in tonus between the two forearm clench muscles?

☒ Yes \_\_\_\_ No

Would you expect to see a difference? Does subject's sex influence your expectations? Explain.

從實驗數據可觀察出慣用手(右手)的力量明顯大於非慣用手，由於慣用手較常受到刺激，藉由去極化和再極化產生動作電位傳遞之後，需要更快速地恢復以便進行下一個動作。

此外，和女性所量測到的數據相比，由於男性的肌纖維較多，因此量測到的肌力較大。

- G. Explain the source of signals detected by the EMG electrodes.  
此訊號的產生是由於肌肉收縮時，肌肉兩端的電位差所產生出來的訊號，而造成肌肉纖維群收縮的原因主要是因為動作電位 (Action Potential) 的傳遞。透過細胞膜內外離子的進出使細胞膜產生電位差而產生「去極化」、「再極化」的現象，其中去極化為細胞膜上鈉離子通道打開使鈉離子進入膜內，造成細胞膜電位上升，膜內外電位差由負轉為正，反之則為再極化。  
肌肉收縮時，這些現象都會產生不同的動作電位，所以EMG所量測的就是這些不同頻率的動作電位所組成的肌電訊號。
- 
- H. What does the term “motor unit recruitment” mean?  
motor unit recruitment可看作為肌肉逐漸增加運動單元(motor unit)來達到充足力量的過程。一個運動單元(motor unit)是由一個motor neuron和附近受它刺激的肌肉纖維群所組成，而肌肉是由數個motor unit組成，肌肉纖維則分散在motor unit之中，當某個motor unit被刺激後，周遭與其相連的unit也會逐漸被刺激起來，因而使肌肉纖維產生足夠的電位。
- 
- I. Define skeletal muscle tonus.  
skeletal muscle tonus 是指在正常情形下，骨骼肌收縮放鬆時仍然會保持在準備狀態而維持一定的輕微張力，並不會一次達到最大力量限度，一般而言此種現象是由中樞神經所控制，使舒張肌肉得以維持正常的長度和張力並保持靈活性，以便進行快速的反應。
- 
- J. Define electromyography.  
Electromyography is a technique which can detect and record the muscle's electrical activities by measuring the voltage changes of the muscle cells when these cells are electrically or neurologically activated with electrodes.
- 

**End of Lesson 1 Data Report**